




V-HOOK

Vitrification Straws



Transparent Sleeve
3 mm Seal



Hook shaped
Flexible Tip



Flat Shaped Handle
for Writing Details

- Hook shaped tip design
- Secure coupling
- Adequate flat surface for writing
- Hook width ensuring minimum placement of solution
- Flexible tip
- 1-cell MEA 80% expand blasto cysts within 96h
- Endotoxin LAL tested
- Sterility: SAL 10^{-6} - Gamma
- Vitrification of oocytes & embryos

Available in
5 colours

Colour	White	Red	Green	Blue	Yellow
Order No.	VH-W	VH-R	VH-G	VH-B	VH-Y

Qty. / Box :
10 Pieces X 2 Pouches

**Surgimedik
Healthcare**

Mobile : +91-771007 9206
Telephone : +91-22-28609401

Domestic : domsales@surgimedik.com
International : sales@surgimedik.com

www.surgimedik.com