

GUIDEWIRES

Schuller Guidewires

For percutaneous entry into vessel using the seldinger technique or for endoscopic access.



- For easy catheter manipulation and introduction of instruments.
- Catheters and other instruments can easily be introduced over the floppy part connected to the rigid shaft.

QTY/Box: 1

Order No.	Shaft	Length (cm)	Rigid Shaft	Floppy Part	Tip
SHJ 100-30-60	.035"	100	30cm	60 cm	J-curved, 3.0mm



SURGIMEDIK