

RGT - Replacement Gastrostomy Tube

Used for replacement gastrostomy placement to provide enteral access for the delivery of nutrition, medication, gastric drainage or decompression.



- Medical grade silicone with soft dome
- Ease of replacement
- Easy traction removal

QTY /Box: 1

Order No.	Description
RGT 16	Replacement Gastrostomy Tube Kit 16 Fr
RGT 20	Replacement Gastrostomy Tube Kit 20 Fr
RGT 24	Replacement Gastrostomy Tube Kit 24 Fr

Each kit includes:

- 1 ea Silicone Gastrostomy Tube with External Bolster
- 1 ea Obturator
- 1 ea Dual Port Feeding Adapter