

Speedo Biopsy Needles

Used for soft tissue biopsy.



- User friendly, lightweight, compact and provided with a strong mechanism : It is therefore more reliable in cutting the most compact tissue
- Provided with 2 triggers: One in the front and one in the back to facilitates operations in every position
- Safety : The safety switch is automatically activated when loading the gun
- Centimeter marking facilitates precise depth placement
- The tip is treated to make it echogenic and facilitate its positioning under ultrasound guidance

QTY. / Box : 5

Order No.	Length (cm)	Gauge	Description
SP16/10	10	16	Speedo, 16 G, 10 cm
SP16/16	16		Speedo, 16 G, 16 cm
SP16/20	20		Speedo, 16 G, 20 cm
SP16/25	25		Speedo, 16 G, 25 cm
SP18/10	10	18	Speedo, 18 G, 10 cm
SP18/16	16		Speedo, 18 G, 16 cm
SP18/20	20		Speedo, 18 G, 20 cm
SP18/25	25		Speedo, 18 G, 25 cm

Domestic:

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.

Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

International:

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com