

### Tumor Stents

The Tumor Stent is used for temporary as well as for palliative treatments of internal urinary drainage when tumors are present in ureteral tract.




- Biocompatible polyurethane with low encrustation tendency and superior radiopacity for enhanced visualization
- Soft and flexible for patient comfort allowing prolonged indwelling period
- Full Shaft with extended diameter

QTY. / Box : 5

Fr	Length (cm)	GW	Order No.
6	24 26 28	0.038"	TS 6/24 TS 6/26 TS 6/28
7	24 26 28	0.038"	TS 7/24 TS 7/26 TS 7/28

### Super Strong Tumor Stents

The Super Strong Tumor Stent is used for temporary as well as for palliative treatments of internal urinary drainage when tumors are present in ureteral tract.



- Soft Kidney & Bladder loops for patient comfort
- 0.035" guidewire compatibility
- Rigid Shaft for maximum compression resistance
- Drainage holes only in the loops
- Radiopaque marker band

QTY. / Box : 5

Fr	Length (cm)	GW	Open Tip
6	24 26 28 30	0.035"	TS 6/24-PA TS 6/26-PA TS 6/28-PA TS 6/30-PA

\* **Combinations** : (With Guidewire : Add Suffix - G)

\* **Each stent is supplied with pusher**

**Domestic:**

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.

Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

**International:**

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com