

Pull-it Stent Remover

Used for Ureteral Stent Removal.




- Effective for various stent materials
- Sharp prongs - firm grip on stent
- Flexible Shaft for easier manipulation

QTY. / Box : 1

Fr	Length (cm)	Order No.
3	60	SR 3/60
4		SR 4/60
5		SR 5/60

Stent Remover Snare

Used for retrieve stents from the urinary tract.



- Nitinol wire is flexible for easy capture
- Nitinol shaft with polymide sheath for trackability

QTY. / Box : 1

Fr	Length (cm)	Order No. (with polymide)
2.4	115	SRS 2.4/115-PA

Domestic:

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.
Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

International:

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com