

## Safety Wire Introducer

Used for Placement of safety wire into the pelvis.



- The outer sheath allows passage of second "Safety" wire during Endourological procedure

QTY. / Box : 10

Fr	Length (cm)	Order No.
10/6	20	SWI

## Nephro Track Nephrostomy Track Balloon Dilator

The Nephrotrack Balloons are recommended for dilating the nephrostomy track.



- Accepts 0.038" Guidewire
- 12 atm rated burst pressure
- 10 mm (30 Fr) X 12 cm balloon
- Supplied with or without Amplatz Sheath
- Option of set with Inflation device 20 cc

QTY. / Box : 1

Balloon Inflated OD	Balloon Length	Shaft Fr	Usable Length	Description	Order No.	Order No. (With Inflation Device)
10 mm (30 Fr)	12	7	55	Without Amplatz	NT 10/12	NT 10/12-IDA
				With Amplatz	NTAM 10/12	NTAM 10/12-IDA

### Domestic:

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.

Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

### International:

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com