


Ultima Percutaneous Nephrostomy Drainage catheter



- A convenient complete set for percutaneous placement of drainage catheter - ultima pigtail design
- Material provides superior biocompatibility and drainage
- Highly radiopaque for precise placement
- Locking loop design available for better anchoring
- Firm polyurethane

QTY. / Box : 10

Fr	Order No. 30 cm	Order No. 30 cm (with UBC)
8	U-PCN-P8/30	U-PCN-P8/30-U
9	U-PCN-P9/30	U-PCN-P9/30-U
10	U-PCN-P10/30	U-PCN-P10/30-U
12	U-PCN-P12/30	U-PCN-P12/30-U
14	U-PCN-P14/30	U-PCN-P14/30-U

Combinations : (With lunderquist Guidewire instead of PTFE Guidewire : Add Suffix - L), (With Hydrophilic Coating : Add Suffix : H)

Ultima Direct Puncture Drainage Set

Used for percutaneous nephrostomy following the "Direct-Puncture" technique particularly when there is a dilated urine collection system a hydronephrosis.



- The Catheter is back-loaded on a two-part Trocar Needle and straightened and then introduced through a skin incision directly
- The procedure does away with the use of guidewire and fascial dilators as there is no dilatation performed
- The set may also be used for cyst and abscess drainage
- with Stop Cock
- Contact us for other sizes

QTY. / Box : 1

Fr	Order No. 30 cm
8	U-DPD-8/30
9	U-DPD-9/30
10	U-DPD-10/30
12	U-DPD-12/30
14	U-DPD-14/30

Combinations : (With Locking Loop System : Add Suffix : LS), (With Hydrophilic Coating : Add Suffix : H)

Domestic:

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.

Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

International:

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com