

## Metal Dilator

Used dilation of penile ureteral meatus.



- Available in adult and pediatric version.

QTY. / Box : 10

Order no.	Description
MDA	Metal Dilator Adult
MDP	Metal Dilator Pediatric

## Disposable Inflation Gauge

Used with balloon dilation catheters to monitor inflation pressure.



- Pressure upto 300 psi
- Cost effective
- Supplied with stop cock to maintain pressure

QTY. / Box : 1

Order no.	Description
PG-V3	Inflation Gauge

## U-Flex Ureteral Access Sheath

Used to establish continuous working channel thereby reduces ureteral trauma during multiple instrument exchange.



- Kink resistant sheath with braided design
- Radiopaque marker bands
- Long Tapered dilator for gradual dilation and less trauma
- Hydrophilic for easier insertion
- Protects delicate instruments & flexible endoscope from damage

QTY. / Box : 1

OD	ID	Length (cm)	Order No.
14 Fr	12 Fr	24	UF 14/12-24
		38	UF 14/12-38
		45	UF 14/12-45
		54	UF 14/12-54
11.5 Fr	9.5 Fr	35	UF 11.5/9.5-35
		45	UF 11.5/9.5-45
		55	UF 11.5/9.5-55

**Domestic:**

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.

Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

**International:**

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com