


## Initial Puncture Needle - 3 Part Bevel Tip

Used for percutaneous nephrostomy and all other percutaneous application where a 0.035" guidewire is to be used.



- Stainless Steel construction
- 3 Part to allow safe puncture
- Bevel Tip

QTY. / Box : 5

Description	Order No.
Initial Puncture Needle, 18 G, 20 cm, 3 Part Bevel	IPN-18-20-3

## Initial Puncture Needle - 2 Part Trocar Tip

Used for percutaneous nephrostomy and all other percutaneous application where a 0.035" guidewire is to be used.



- Stainless Steel construction
- Trocar Tip

QTY. / Box : 5

Description	Order No.
Initial Puncture Needle, 18 G, 20 cm, 2 Part Trocar	IPN-18-20-2T
Initial Puncture Needle, 18 G, 30 cm, 2 Part Trocar	IPN-18-30-2T

**Domestic:**

Surgimedik Healthcare India Pvt. Ltd., Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.  
Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404 domsales@surgimedik.com

**International:**

Surgimedik Healthcare, 223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA. sales@surgimedik.com