


# Flexi Cath C

2



## Curved Dual Lumen Catheter

- Soft tip for easier insertion
- Bodysoft polyurethane material
- One step “Stepper Dilator” included
- Available in catheter & kit version
- Seldinger technique insertion
- Curved Extention for comfort and ease of maintenance

Colour coded luers



Soft Tip



### Catheter Contents :

1 ea Flexi-Cath-C, 1 ea Stepper Dilator 12 F, 2 ea Injection Caps.

### Kit Contents :

1 ea Flexi-Cath-C, 1 ea Stepper Dilator 12 F, 2 ea Injection Caps.

1 ea Introducer Needle 18/7, 1 ea Guidewire JFS 80-38, 1 ea 5ml Syringe, 1 ea Scalpel.

Order No.	Description	Qty. / Box
DLCC 12/13	Flexi-Cath-C Curved Dual Lumen Catheter 12 F, 13 cm	1
DLCC 12/16	Flexi-Cath-C Curved Dual Lumen Catheter 12 F, 16 cm	1
DLCC-K 12/13	Flexi-Cath-C Curved Dual Lumen Catheter Kit 12 F, 13 cm	1
DLCC-K 12/16	Flexi-Cath-C Curved Dual Lumen Catheter Kit 12 F, 16 cm	1

### **Domestic:**

Surgimedik Healthcare India Pvt. Ltd.  
Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.  
Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404  
domsales@surgimedik.com

### **International:**

Surgimedik Healthcare,  
223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA.  
sales@surgimedik.com



**SURGIMEDIK**

Patient's life is our life