

PD Cath

5

Stylet Catheter for Peritoneal Dialysis

Precision Hole Punching:

Allows sufficient flow for quick administration and drainages

Elbow Connector:

Prevents catheter kinking to ease application

Catheter Marking:

Shows direction of distal end for facilitate insertion

Roller Clamp:

Facilitate tubing open / close procedure for ease procedure

Rubber Connector:
Self Sealing to allows easy additional medication

Tapered Catheter Tip:
Provides ease of placement to patient comfort

*Available in Adult & Pediatric versions.


Thumbring for Better Grip




Trocar Tip for easier inserstion



Curved Catheter Connector



Movable Ball



Flow Regulator



Catheter Contents :

1 ea PD Catheter, 1 ea Catheter Connector, 1 ea Scalpel.

Order No.	Description	Qty. / Box
PDCA	PD Catheter Adult 30 cm	1
PDCP	PD Catheter Pediatric 20 cm	1
PDX	PD Transfer Set	1

Domestic:
Surgimedik Healthcare India Pvt. Ltd.
Surgimedik House, 53 CD, Kandivali Industrial Estate, Charkop, Kandivali (West), Mumbai - 400067, INDIA.
Tel.: +91 22-2860 9401 Fax : +91 22-2860 9404
domsales@surgimedik.com

International:
Surgimedik Healthcare,
223, Surat Special Economic Zone, Sachin, Surat - 394230, Gujarat, INDIA.
sales@surgimedik.com



SURGIMEDIK

Patient's life is our life