


EFFICIENT

Intra Uterine Insemination Catheter



Used for introduction of washed spermatozoa into the uterus.

- Soft Rounded Tip with Hole for atraumatic placement
- Adjustable Clear Silicone Stopper to visualise insertion depth of catheter
- Inner Flexible Sheath (Soft Tubing) for patient comfort & smooth cervical passage
- Outer Sheath (Hard Tubing) for easy trackability & ease of insertion
- Non-toxic medical grade tubing
- LAL tested
- Without centimeter marking



Order No.	Qty./Box
E-IUI	50