



R+ Multiloop Ureteral Stents



The R+ Multiloop Ureteral Stent is intended to temporarily facilitate drainage of urine from kidney to bladder by trained physician.


- Multiloop designs accommodates ureteral lengths from 22 cm to 30 cm
- With enhanced radiopacity
- Contact us for other sizes



Pelvis Side
Open Tip



Pelvis Side
Closed Tip



Centimeter Marking
on Stents

Stent + Pusher

Qty. / Box : 10 or 25

Open Tip (Bladder End Multiloop) Order No.	Closed Tip (Bladder End Multiloop) Order No.	Open Tip (Both End Multiloop) Order No.	Closed Tip (Both End Multiloop) Order No.	Fr	Length (cm)	Compatible Guidewire
SOT 4.7 / ML - R+	SS 4.7 / ML - R+	SOT 4.7 / ML - ML - R+	SS 4.7 / ML - ML - R+	4.7	22-30	0.035"
SOT 6 / ML - R+	SS 6 / ML - R+	SOT 6 / ML - ML - R+	SS 6 / ML - ML - R+	6	22-30	0.038"
SOT 7 / ML - R+	SS 7 / ML - R+	SOT 7 / ML - ML - R+	SS 7 / ML - ML - R+	7	22-30	0.038"

** Set available with guidewire, suture, clamp etc.
Hydrophilic coating on stent as per requirement.