



EAGLE-FL

Single Lumen Ovum Pickup Needles with Flushing Line



Used for laparoscopic or ultrasound guided transvaginal harvesting of oocyte from follicles along with flushing of follicles.



- Highly echogenic design
- Silicone bung for tight fit
- Non-coring needle
- MEA & LAL tested
- With luer lock connector

- Regulation of flow
- Alternate aspiration & flushing made possible through single channel
- Multiple operation functionality by using stop cock valve
- Available in 17 G

Order No.	Gauge	Length	Qty./Box
E-OPN-SL17/35-FL	17	35 cm	10