



V-STRIP

Vitrification
Straws



Transparent Sleeve
3 mm Seal








Flexible Transparent
Strip with tip marking



Flat Shaped Handle
for writing details

- High survival and development rates.
- Secure coupling.
- Marking on tip extreme for easy visibility under microscope.
- Adequate flat surface for writing.
- Tip width ensuring placement of minimum solution
- 1-cell MEA \geq 80% expand blastocysts within 96h
- Endotoxin LAL tested
- Sterility: SAL 10^{-6} - Gamma
- Vitrification of oocytes & embryos in all stages of development
- Available in 5 colours

Available in
5 colours

					
Colour	White	Red	Ivory	Blue	Yellow
Order No.	VS-W	VS-R	VS-I	VS-B	VS-Y

Qty. / Box :
10 Pieces X 2 Pouches

**Surgimedik
Healthcare**

Mobile : +91-771007 9206
Telephone : +91-22-28609401

Domestic : domsales@surgimedik.com
International : sales@surgimedik.com

www.surgimedik.com