


V-LOCK Vitrification Straws



Semi Transparent Tip with marking




Flat Shaped Handle for Writing Details








Semi Transparent Cap

- Adequate safety
- Secure air tight coupling.
- Semi transparent body and cap.
- No extra accessories required.
- Arrow on top surface ensuring orientation during use.
- Shallow well on tip distal end for exact placement of embryo.
- Ease in visibility under microscope
- 1-cell MEA \geq 80% expand blastocysts within 96h
- Endotoxin LAL tested
- Sterility: SAL 10^{-6} - Gamma
- Vitrification of in 1-cell stage of embryos

Available in **5 colours**



					
Colour	Pink	Violet	Green	Blue	Yellow
Order No.	VL-P	VL-V	VL-G	VL-B	VL-Y

Qty. / Box :
10 Pieces X 2 Pouches